

## SeromaCath® WOUND DRAINAGE SYSTEM

*specifically designed to treat postoperative seromas*

Greer Medical's patented SeromaCath® Wound Drainage System is the first and only wound drainage system specifically designed and clinically proven to effectively treat postoperative seromas following: mastectomy, axillary dissection, groin dissection, breast reconstruction and augmentation, abdominoplasty, liposuction, latissimus dorsi flap, TRAM flap, and other procedures.

The newest in our family of catheters is the BLUNT SeromaCath® specifically designed to treat seromas surrounding breast implants without the risk of perforation.

U.S. Patent No. 4,341,212

# SeromaCath®

WOUND DRAINAGE SYSTEM

**Greer**  
MEDICAL, INC.

CE 0336 ISO 9002/EN 46002

[www.greer-medical.com](http://www.greer-medical.com)

# SeromaCath®

WOUND DRAINAGE SYSTEM

*The most effective method for treating a variety of postoperative seromas following:*

- *mastectomy*
- *axillary dissection*
- *groin dissection*
- *breast reconstruction and augmentation*
- *abdominoplasty*
- *liposuction*
- *latissimus dorsi flap*
- *TRAM flap*
- *and other procedures*

**ORIGINAL  
cat. #1002**

14G 2" Teflon® catheter/stylet needle/plug with six side ports and a connecting drainage tube for a bulb suction reservoir.

1 per package/5 packages per box

**LONGER  
cat. #1035**

14G 3.5" Teflon® catheter/stylet needle/plug with six side ports and a connecting drainage tube for a bulb suction reservoir.

1 per package/5 packages per box

**BLUNT  
cat. #1275**

14G 2" Teflon® catheter/stylet needle/plug with two side ports, a 2.75" blunt metal catheter with two side ports and a connecting drainage tube for a bulb suction reservoir.

1 per package/5 packages per box

**RESERVOIR  
cat. #2001**

100cc bulb suction reservoir with anti-reflux valve and clip.

1 per package

**To place your order call 1-800-424-2155**

**Greer**


MEDICAL, INC.

314 East Carrillo Street  
Santa Barbara, California 93101 USA  
Telephone: 805-962-5883  
Fax: (805) 962-4615  
E-mail: greer-medical@msn.com


Authorized Representative in EU:  
EMDAR BV  
IJsselburcht 3  
Postbus 5486  
6802 EI Arnhem  
The Netherlands



**ORIGINAL (1002) SeromaCath®**  
for treating a variety of postoperative seromas.



**LONGER (1035) SeromaCath®**  
for treating postoperative seromas beneath a thick panniculus.



**BLUNT (1275) SeromaCath®**  
for treating postoperative seromas surrounding breast implants.

CE 0336 ISO 9002 / EN 46002 PN: 1002, 1035, 1275 Rev.8/02

U S A • E U R O P E • A S I A